



105mm

40mm

210mm



OPERATION MANUAL

NON-CONTACT DIGITAL INFRARED THERMOMETER



Manufacturer:
Hunan Zhongliankang Medical Technology Co., Ltd.
Factory address:
Second floor, Building 10, Innovation and Entrepreneurship Center, No. 9 Chuanqi West Road, Jiuhua Economic Development Zone, Xiangtan City, Hunan



MEASUREMENT INSTRUCTION

- Align the thermometer probe in the middle of the forehead and keep it vertical with a distance of 3-5 cm. Press the measurement button, and the temperature will be read and displayed after 1 second with sound of 'Tick'. And then, the measurement is completed.
- When the great change of the ambient temperature affects the forehead temperature, please measure the earlobe at its back.



MAINTENANCE AND PRECAUTIONS

- Please keep the inner cavity of the sensor and probe clean, otherwise it will affect the measurement accuracy.**
- Surface cleaning: Wipe surface dirty with a clean soft cloth or cotton swab with a little medical alcohol or water.
 - Sensor and inner cavity cleaning: use a clean soft cloth or cotton swab to stick a little medical alcohol gently wipe the inner cavity of the probe or the top of the sensor, use until the alcohol is completely evaporated.
- Please read the instructions carefully before use.
 - Do not immerse the thermometer in any liquid, and do not leave it under overheat or low temperature environment for a long time.
 - Collision prohibited. Dropping and mixing with sharp objects are prohibited, disassemble by yourself is prohibited.
 - The thermometer should not be used in the sun or water.
 - Do not use in strong electromagnetic interference environment.
 - Keep the thermometer out of the reach of children.
 - It is recommended to practice several times to become familiar with the measuring method, and try not to change the product factory settings.
 - Measurement results are not a substitute for physician diagnosis.
 - No special maintenance is required during use, please contact the seller or manufacturer if there is any malfunction appears.

REMARKS

- Before and after use, keep the inner part of the sensor and probe clean.
- Please use the thermometer in a stable temperature environment. When the ambient temperature fluctuates greatly (such as the temperature difference from indoor to outdoor), please wait for 10 seconds before measuring.
- Do not start measuring body temperature immediately after measuring the objects with extremely high or low temperature. Please leave it for 10 minutes before measuring the body temperature.
- When the measured people or object comes from a place with a large temperature difference from the measurement environment, they should stay in the measurement environment for at least 5 minutes.
- Try not to measure when there's blowing, sweating, water, makeup or etc. on the forehead. Don't measure your body temperature in 30 minutes after exercise, bath or meal.
- Press and hold the SET button in power off state, and enter the setting mode after the 'beep' sound. Press the SET button to select F3, and use the MEMO button switches between °C mode and °F mode. After setting, press the POW button to power off.

Product model	ZLK-IRT101
Display resolution	(32.2°F)/(0.1°C)
Body temperature mode	(±32.5°F)/(±0.3°C)
Measuring range	(89.6°F~95°F, ±32.5°F)/(32°C~35°C, ±0.3°C) (95°F~107.6°F, ±32.3°F)/(35°C~42°C, ±0.3°C) (95°F~109.4°F, ±32.5°F)/(35°C~43°C, ±0.3°C)
Measuring distance	1~10cm
Operating ambient temperature	(50°F~104°F)/(10°C~40°C)
Operating ambient humidity	≤80%RH
Storage temperature range	(-13°F~131°F)/(-25°C~55°C)
Power	DC3V
Battery	2 AAA batteries
IP rating	IP 20
Dimension	116x55x170mm(LxWxH)
Weight	125g (without battery)

WARRANTY CARD

Purchaser Name	
Telephone Number	
Product Model	
Company Name	
Date of Purchase	
Distributor's Seal	

CUSTOMER FORM

Purchaser Name	
Telephone Number	
Product Model	
Company Name	
Date of Purchase	
Distributor's Seal	

DISTRIBUTOR FORM